PL - Derma roller

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EN - Derma roller

Usage: Disinfect the roller before the procedure. Thoroughly cleanse and dry the skin before use. Then roll over the skin 10 times in each direction (up, down, sideways, and diagonally). One session should last 3–10 minutes depending on the intensity. To enhance results, a nourishing serum can be applied.

Warnings: The roller is for personal use only. After treatment, redness and small scabs may appear. Avoid sun exposure, tanning beds, sauna, swimming pool, and alcohol before and after the procedure. Not recommended for pregnant or breastfeeding women, those planning pregnancy, tanning bed users, people exposed to sunlight or using self-tanners. Do not use directly on tattoos, dark brown or black spots. Do not use around nipples, genitals, or anus. The number of sessions depends on individual skin condition. Do not use on diseased or damaged skin. After each use, disinfect the roller. Do not place in boiling water. This is not a toy. Use only as intended. Keep out of reach of children and people with physical or mental disabilities.

Contraindications: Contagious skin diseases, inflammation, blood clotting disorders, diabetes, eczema, psoriasis, jaundice, leukemia, herpes, couperose skin, very dry skin, allergies, cancer, fillers or botox.

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