Producent/Manufacturer/Hersteller/Valmistaja/Tootja/Ražotājs/Gamintojas/Výrobce/Výrobca/Gyártó/Producător/Производитель
n VITOLOG. PL sp. z o.o.
ul. Międzyleska 6
50-514 Wrocław
Nip:8992838703,
Kontakt: +48 71 727 62 28
kontakt@vitolog.pl

Wyprodukowano w Chinach/Made in China/Hergestellt in China/Valmistettu Kiinassa/Valmistatud Hiinas/Ražots
Ķīnā/Pagaminta Kinijoje/Vyrobeno v Číne/Vyrobén v Číne/Kínában készült/Fabricat în China/Произведено в Китае

PL - Derma roller 3w1

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EN - Derma Roller 3-in-1

Product Features: Derma roller; material: stainless steel, ABS; 3 interchangeable heads: 300c, 720c, 1200c; device for quick and easy micro-needle mesotherapy; indications: wrinkle removal, cellulite, discoloration, poor skin condition, sagging skin, hair loss.

How to Use: Select the appropriate head and attach it to the handle. Insert it into the handle until you hear a click. To remove the head, press the button on the handle and gently pull it out. Before the procedure, disinfect the roller. Thoroughly cleanse and dry the skin. Roll the skin 10 times in each direction (up, down, sideways, diagonally). One treatment should last 3-10 minutes, depending on the intensity of rolling. To enhance the effect, you can apply a nourishing product.

Warnings: One roller should only be used by one person. Redness and small "scabs" may occur after the procedure. Avoid tanning, solariums, saunas, swimming pools, and alcohol consumption before and after the procedure. The product is not recommended for pregnant or breastfeeding women, those planning pregnancy, solarium users, those who sunbathe, or those using self-tanners. Do not use the product directly on tattoos, dark brown, or black spots. Use around nipples, genitals, or anus is prohibited. The number of treatments depends on the individual skin condition. Do not use on diseased skin. Avoid using IRRITATING cosmetics before the procedure. After each treatment, rinse the roller with disinfectant. Do not place the product in boiling water. The product is not a toy. Use the product as intended. Keep the device away from children.

Contraindications: Infectious skin diseases, inflammation, blood clotting disorders, diabetes, eczema, psoriasis, jaundice, leukemia, herpes, vascular skin, dry skin, allergies, cancer, fillers, and Botox.

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